

meditest

Reagents Catalog



2026

Contents

Meditest, 2026, Reagents Catalog

01

Introduction

02

Biochemistry
All List

03

Biochemistry
Detailed List

14

Hematology
All List

16

Coagulation
All List

18

Serology
All List

20

Blood Grouping
All List

INTRODUCTION

Meditest, 2026, Biochemistry Reagent Product Catalog



Medios is one of the companies specialized in R&D and production of IVD products in Türkiye.

GET TO KNOW US

Medios was founded in 2017. Since its establishment, it has continued its R&D investments without slowing down. We have created a state-of-the-art R&D and production center that produces Microbiology, Biochemistry, Hemogram, Coagulation and External Quality Control materials.

1000 m2 closed area

500 m2 production area

200 m2 R&D area

150+ products

19+ staff



Biochemistry

Reagent Products

- | | | | |
|------------------------|------------------|---------------------|-----------------|
| • Acid Phosphatase | • Cholesterol | • Glucose | • Total Protein |
| • Albumin | • Cholinesterase | • HbA1c | • Triglycerides |
| • Alkaline Phosphatase | • CK-MB | • HDL Cholesterol | • Urea |
| • Ammonia | • CK-NAC | • Iron | • Uric Acid |
| • Amylase | • Complement 3 | • UIBC | • Zinc |
| • Anti-Streptolysin O | • Complement 4 | • LDH | • Vitamin D |
| • AST | • Copper | • LDL Cholesterol | |
| • Bilirubin Total | • Creatinine | • Magnesium | |
| • Bilirubin Direct | • CRP | • Phosphorus | |
| • Calcium | • Gamma GT | • Rheumatoid Factor | |

Calibrators

- Diacheck Biochemistry Calibrator
- Diacheck HDL/LDL Calibrator
- Diacheck HbA1c Calibrator
- Diacheck VitD Calibrator
- Diacheck ASO Calibrator
- Diacheck CRP Calibrator
- Diacheck RF Calibrator
- Diacheck Zinc Calibrator
- Diacheck Copper Calibrator

Controls

- Diacheck Biochemistry Control Set
- Diacheck HbA1c Control Set
- Diacheck VitD Control Set
- Diacheck Protein Control Set
- Diacheck Zinc Control Set
- Diacheck Copper Control Set

Sarflar


- Basic Wash Solution
- Acid Wash Solution
- Hitergent 59mL

Dedicated Reagents for Analyser

- Abbot
 - Architect c4000/c8000
- Beckman Coulter
 - AU 400/480/640/680/2700
- Siemens
 - Advia 1200/1800
- Mindray
 - BS 380/480/800
- Erba
 - XL 600/640/800/1000



Ready Barcode



Original Bottle



Ready Application

Acid Phosphatase Liquid ready-to-use reagents Ref: 24ACP01

Measuring Range : 0,8-150 U/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

Albumin BCG Liquid ready-to-use reagents Ref: 24ALB01

Measuring Range : 0,01-6 g/dL

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

ALP p-Nitrophenylphosphate Liquid ready-to-use reagents Ref: 24ALP01

Measuring Range : 7-3000 U/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

ALT NADH. Kinetic UV. IFCC Liquid ready-to-use reagents Ref: 24ALT01

Measuring Range : 0-400 U/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

Ammonia (Enzymatic UV) Liquid ready-to-use reagents Ref: 24AMN01

Measuring Range : 26,2-600 umol/L

Calibrator: Diacheck Ammonia Calibrator

Control : Diacheck Ammonia Control Set

Stability : 28 day on board the analyser at 2-8 °C

Amylase CNPG3. Kinetic

Liquid ready-to-use reagents **Ref: 24AMY01**

Measuring Range :	0,91-2000 U/L
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Anti-Streptolysin O Turbidimetric

Liquid ready-to-use reagents **Ref: 24ASO01**

Measuring Range :	3-800 IU/mL
Calibrator:	Diacheck ASO Calibrator
Control :	Diacheck Protein Control Set
Stability :	28 day on board the analyser at 2-8 °C

AST NADH. Kinetic UV. IFCC

Liquid ready-to-use reagents **Ref: 24AST01**

Measuring Range :	3-800 U/L
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Bilirubin Direct DMSO

Liquid ready-to-use reagents **Ref: 24BILD01**

Measuring Range :	0,01-13 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Bilirubin Total DMSO

Liquid ready-to-use reagents **Ref: 24BILT01**

Measuring Range :	0,03-25 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Calcium Arsenazo III Liquid ready-to-use reagents Ref: 24CA01

Measuring Range :	0,024-20 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Cholesterol CHOD-POD Liquid ready-to-use reagents Ref: 24CHOL01

Measuring Range :	0,3-1000 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Cholinesterase Liquid ready-to-use reagents Ref: 24CHE01

Measuring Range :	432-2500 U/L
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

CK-MB Immunoinhibition Liquid ready-to-use reagents Ref: 24CKMB01

Measuring Range :	10-2000 U/L
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

CK-NAC IFCC Liquid ready-to-use reagents Ref: 24CKNAC01

Measuring Range :	10-2000 U/L
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Complement C3 Turbidimetry

Liquid ready-to-use reagents

Ref: 24C301

Measuring Range :	13-502 mg/dL
Calibrator:	Diacheck Protein Calibrator
Control :	Diacheck Protein Control Set
Stability :	28 day on board the analyser at 2-8 °C

Complement C4 Turbidimetry

Liquid ready-to-use reagents

Ref: 24C401

Measuring Range :	2,90-152 mg/dL
Calibrator:	Diacheck Protein Calibrator
Control :	Diacheck Protein Control Set
Stability :	28 day on board the analyser at 2-8 °C

Copper

Liquid ready-to-use reagents

Ref: 24COP01

Measuring Range :	4-500 ug/dL
Calibrator:	Diacheck Copper Calibrator
Control :	Diacheck Copper Control Set
Stability :	28 day on board the analyser at 2-8 °C

Creatinin Jaffe Modified Method

Liquid ready-to-use reagents

Ref: 24CREA01

Measuring Range :	0,20-20 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

CRP Turbidimetry

Liquid ready-to-use reagents

Ref: 24CRP01

Measuring Range :	0.5-175 mg/L
Calibrator:	Diacheck CRP Calibrator
Control :	Diacheck Protein Control Set
Stability :	28 day on board the analyser at 2-8 °C

GGT Carboxy substrate. Kinetic

Liquid ready-to-use reagents

Ref: 24GGT01

Measuring Range : 2-300 U/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

Glucose GOD-POD

Liquid ready-to-use reagents

Ref: 24GLU01

Measuring Range : 1-500 mg/dL

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

HbA1c Immunoturbidimetric

Liquid ready-to-use reagents

Ref: 24A1c01

Measuring Range : 4-16%

Calibrator: Diacheck A1c Calibrator

Control : Diacheck A1c Control Set

Stability : 28 day on board the analyser at 2-8 °C

HDL Direct

Liquid ready-to-use reagents

Ref: 24HDL01

Measuring Range : 5-200 mg/dL

Calibrator: Diacheck Lipid Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

IRON Ferrozine

Liquid ready-to-use reagents

Ref: 24FE01

Measuring Range : 25-1000 ug/dL

Calibrator: Diacheck Iron Standart

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

UIBC Ferrozine Liquid ready-to-use reagents Ref: 24UIBC01

Measuring Range : 19-700 ug/dL

Calibrator: Diacheck Iron Standart

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

LDH Pyruvate. Kinetic UV. DGKC Liquid ready-to-use reagents Ref: 24LDH01

Measuring Range : 31-4000 U/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

LDL Direct Liquid ready-to-use reagents Ref: 24LDL01

Measuring Range : 4-500 mg/dL

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

Magnesium Xylidyl Blue. Colorimetric Liquid ready-to-use reagents Ref: 24MG01

Measuring Range : 0,036-3,27 mmol/L

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

Phosphorus Liquid ready-to-use reagents Ref: 24PHOS01

Measuring Range : 0,4-20 mg/dL

Calibrator: Diacheck Biochemistry Calibrator

Control : Diacheck Biochemistry Control Set

Stability : 28 day on board the analyser at 2-8 °C

RF Turbidimetry

Liquid ready-to-use reagents

Ref: 24RF01

Measuring Range :	2-160 IU/ml
Calibrator:	Diacheck RF Standart
Control :	Diacheck Protein Control Set
Stability :	28 day on board the analyser at 2-8 °C

Total Protein Biuret

Liquid ready-to-use reagents

Ref: 24TP01

Measuring Range :	0,10-12 g/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Triglycerides GPO-POD

Liquid ready-to-use reagents

Ref: 24TRIG01

Measuring Range :	1-1600 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Urea Urease-GLDH. Kinetic

Liquid ready-to-use reagents

Ref: 24UREA01

Measuring Range :	1.24-530 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

Uric Acid Uricase-POD. Liquid

Liquid ready-to-use reagents

Ref: 24UA01

Measuring Range :	0,02-40 mg/dL
Calibrator:	Diacheck Biochemistry Calibrator
Control :	Diacheck Biochemistry Control Set
Stability :	28 day on board the analyser at 2-8 °C

ZINC Colorometric Liquid ready-to-use reagents Ref: 24ZN01

Measuring Range : 11.3 - 159 $\mu\text{mol/l}$

Calibrator: Diacheck ZINC Calibrator

Control : Diacheck ZINC Control Set

Stability : 28 day on board the analyser at 2-8 °C

Vitamin D ImmunoTurbidimetric Liquid ready-to-use reagents Ref: 24VD01

Measuring Range : 4-150 ng/mL

Calibrator: Diacheck Vitamin D Calibrator

Control : Diacheck Vitamin D Control Set

Stability : 28 day on board the analyser at 2-8 °C

Roche Hitergent ECO-D Solution Liquid ready-to-use reagents Ref: 24ECO01

Stability : 28 day on board the analyser Packing : 6x10 mL

Roche Basic Wash Solution Liquid ready-to-use reagents Ref: 24NAH01

Stability : 28 day on board the analyser Packing : 2x2 L

Roche Acid Wash Solution Liquid ready-to-use reagents Ref: 24CAS01


Stability : 28 day on board the analyser Packing : 2x2 L

Dedicated Original Bottle and Ready Barcode

Abbott Series



Ready Barcode



Original Bottle



Ready Application

Beckman Coulter Series



Ready Barcode



Original Bottle




Ready Application


Siemens Series



Ready Barcode



Original Bottle




Ready Application

Dedicated Original Bottle and Ready Barcode

meditest

Roche Series



Clean Solutions
and Reagent

Ready Application

Mindray Series



Ready Barcode



Original Bottle




Ready Application

Dirui Series



Ready Barcode



Original Bottle



Ready Application

Hematology

18 Parameter Reagent Products

- Diluent Solution
- Lyse Solution
- Rinse Solution
- Probe Cleanser
- 3 DIFF Control Set (L,N,H)



22 Parameter Reagent Products

- Diluent Solution
- Baso Lyse Solution
- RBC Lyse Solution
- Rinse Solution
- Probe Cleanser
- 5 DIFF Control Set (L,N,H)



Applicable Analysers

- Mindray
 - BC3000 Plus
 - BC30s
 - BC6000 Series
 - BC700 Series
- Erma
- Dymind
 - UN73
- Sysmex
 - Kx21



Capacity to produce
50 million tests per month

Coagulation

Reagent Products

- Recombi Pt Reagent
- PT Reagent
- Aptt Reagent
- CaCl₂ Reagent
- Fibrinogen Reagent

Applicable Analysers

- Roche t411
- Zonci Series
- Succeder Series
- Werfen Series
- Rayto Series
- Alaris
- And other analysers

Prothrombin Time Reagent

Cat No	Product	Packing	ISI	Source	Open Stability
24rPT01	Recombi PT-s	10x4 mL	1.0-1.3	Human	10 Day
24rPT02	Recombi PT-s	10x20 mL	1.0-1.3	Human	10 Day
24PT01	PT-s	10x4 mL	1.0-1.3	Rabbit	7 Day
24PT02	PT-s	10x20 mL	1.0-1.3	Rabbit	7 Day

Active Partial Thromboplastin Time Reagent

Cat No	Product	Packing	Activator	Format	Open Stability
24APTT01	aPTT-s	10x4 mL	Ellagic Acid	Liquid	30 Day
24APTT02	aPTT-s	10x20 mL	Ellagic Acid	Liquid	30 Day
24CC01	Calcium Chloride	10x4 mL	-	Liquid	Expiry Date
24CC01	Calcium Chloride	10x20 mL	-	Liquid	Expiry Date

Fibrinogen



Cat No	Product	Packing	Format	Open Stability
24FIB01	Fibrinogen-s	10x4 mL	Liquid	10 Day
	Imidazol Buffer	10x20 mL	Liquid	7 Day
	Fibrinogen Cal	10x4 mL	Liquid	6 hours

Serology

Reagent Products

- Brucella Rose bengal with control
- Brucella abortus
- Brucella melitensis
- Salmonella typhi spp.
- Salmonella paratyphi spp.
- Aso Latex
 - Aso control
- Crp Latex
 - Crp control
- RF Latex
 - Rf control

Cat No	Product	Packing	Shelf Time
25RBL01	Meditest Rose Bengal 100 Test	1 x 5 mL	18 Month
25RBL02	Meditest Rose Bengal 200 Test	1 x 10 mL	18 Month
25RBL03	Meditest Rose Bengal 100 Test with Control	1 x 5 mL	18 Month
25BAB01	Meditest Brucella abortus	1 x 5 mL	18 Month
25MEL01	Meditest Brucella melintesis	1 x 5 mL	18 Month
25SAT01	Meditest Salmonella typhi spp	1 x 2 mL	18 Month
25SAP01	Meditest Samonella paratyphi spp	1 x 2 mL	18 Month



Blood Grouping

Reagent Products

- Anti A Monoclonal
- Anti B Monoclonal
- Anti AB Monoclonal
- Anti D Monoclonal
- Bovine Albumin %22
- Quality Control Set
- Liss
- Bromeline

Cat No	Product	Packing	Shelf Time
25BG01A	Meditest Anti A Monoclonal	1 x 10 mL	24 Month
25BG01B	Meditest Anti B Monoclonal	1 x 10 mL	24 Month
25BG0AB	Meditest Anti AB Monoclonal	1 x 10 mL	24 Month
25BG01D	Meditest Anti D Monoclonal	1 x 10 mL	24 Month
25BGBVN	Meditest Bovine Albumine %22	1 x 10 mL	18 Month
25BG0QC	Meditest Quality Control Set	4 x 1 mL	2 Month
25BGLISS	Meditest LISS Solution	1 x 500 mL	24 Month
24GBBML	Meditest Bromeline	1x5 mL	24 Month



meditest

REAGENT INDUSTRY

“We are always working to make it better”

"Medios aims to create value for its employees, suppliers and customers and to become one of the leading players in the field of diagnostics in the world by transferring this value to future generations."

Contact:

Adress: Guzelevler Mahallesi Cinar Caddesi No: 5
Yuregir/Adana

Tel: 0322 502 0095

Web: www.medios.web.tr - www.medios.com.tr

e-Mail: medios@medios.com.tr