

HEMOGRAM

02

nexus
rc360

- You get 60 tests.
- HD color screen, sensitive touch.
- Detects touch with gloves.
- Cyanide-free HGB measurement technology

**Intelligent measuring
technology combines high
performance and compact
design**

NEXUS rc360
60 TESTS/TIME
<10 µL Sample volume



- 18 parameters
- 10 µL then as sample
- High accuracy in results
- 10.4 inch color screen

Measurement Parameters: WBC, Lym#, Mid#, Gran#, Lym%, Mid%, Gran%, RBC, HGB,HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT, P-LCR, P-LCC

HEMOGRAM

nexus
rc560

- Per hour 60 tests.
- Laser Scattering + Sheat Flow Technology
- Emergency sample position.

Industry's smallest 50 μm
diameter for unparalleled
precision

NEXUS rc560
60
TESTS/HOUR
<10 μL Sample volume



- 25 parameters
- 10 μL then as sample
- >100,000 patient storage
- 10.4 inch color screen

Measurement Parameters: WBC, LYM%, LYM#, NEU%, NEU#, MON%, MON#, EOS%, EOS#, BAS%, BAS#, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, PDW, MPV, PCT, P-LCR, P-LCC